

Nordic guidelines for dose reduction to radiosensitive organs of the patient in conventional radiography and fluoroscopy

Anders Widmark, Anders Beierholm, Carl Bladh, Ritva Bly, Eva G. Friberg, Nelly Pétursdóttir, Petra Tenkanen-Rautakoski, Hanne N. Waltenburg

Statement from The Nordic Radiation Protection co-operation

Nordic Working Group on Medical Applications



Anders Widmark
anders.widmark@dsa.no

Norwegian Radiation
and Nuclear Safety Authority

The Nordic Radiation Protection co-operation for the radiation protection authorities in medical applications

- 40 years celebration last year
- Previous statements <http://nordicxray.gr.is>
 - Concerning the increase use of computed tomography in the Nordic countries
 - Bismuth shielding of patients in CT examinations
 - Position statement on justification of new types of practices involving medical exposure
 - See also *Session S6-P1 History of NGMA*
- The current statement discuss the optimization in radiography for proper use of equipment and radiographic technique in medical applications

Why a statement for basic radiography?



- Maybe education and training for radiographers for CT and MRI is favored in relation to basic radiography?
- Indication of high rates of re-takes in the collection of local DRL

Premises for this statement

- A functional system for QA and QC
- AEC is properly adjusted
- Filtration is optimized for the actual procedure and laboratory



Projections – AP vrs. PA

- HVL in human tissue = ~3 cm in conventional radiography (75 % of the dose in first 6 cm)
- Normally ~1 % of the incident radiation reach detector
- Organs more on the ventral (front) side: eye lens, breast tissue, intestines, uterus... will have a significant dose reduction in PA projections
- ICRP 2012 – Concern about radiosensitivity of the eye lens
- ICRP 2007 – W_T increased for breast tissue from 0.05 to 0.12
- W_T factors are an average over the whole population (all ages and both sex)
- Special precautions for young females and exposure of breast tissue in mammography
- *PA projections - A significant decrease in radiation dose to many sensitive organs*



Collimation of radiation field to ROI

- Avoid irradiating other organs unnecessary
- Reduces the need for shielding of radiosensitive organs
- Reduces the amount of scattered radiation – Scatter increases linear with the irradiated area
- Rectangular collimator in dentistry (requirement in Norwegian regulation from 2020)
- *Strict collimation - Reduced scatter to other organs in patient, decreased staff exposure and increased image quality*



Use of grid

- Increase the image quality, but also the dose $\sim x 3$
- Usually not necessary for children, due to smaller exposed volume
- When planning for paediatric laboratory, look for equipment where the grid can be easily taken away, also for fluoroscopy
- Consider if the grid is necessary when imaging children

Compression



- More common in the past, but got a new renaissance due to more effective compression equipment
- Most suitable for pelvis, lumbar spine and non-acute abdomen
- HVL in tissue ~3 cm
- Most patients can be compressed 7-8 cm in the abdominal area
- Also, avoiding movement unsharpness
- *Consider to implement the use of compression*

Gonad shielding – General comments

- Very important in the 60's and 70's
- Used mainly for reducing the risk for hereditary effects
- ICRP 2007 reduced the risk estimates for hereditary effects by a factor of 6-8 times
- Reflected in the decrease of W_T for gonads, from 0.20 to 0.08
- One have to admit: Shielding of gonads was more important in the past!
- UNSCEAR: «No radiation-induced hereditary diseases have so far been demonstrated in human populations. However, experimental studies in plants and animals have clearly demonstrated that radiation can induce genetic effects; consequently, humans are unlikely to be an exception in this regard.»




Gonad shielding - Males

- Re-cycling of sperm have a cyclus of about 70 days
- Shielding - Option when the radiation field is ≤ 5 cm from testes
- Dedicated shielding equipment (capsules), age and size specific
- May reduce testes dose up to 95 %, if done properly
- If it is a risk for re-takes or interfering with AEC - no shielding
- *Use of gonad shielding on males are effective when used properly*



Gonad shielding - Females

- The position of the ovaries can vary significantly, especially for young females
- Contact shielding will not shield from scatter produced inside the body
- Max. 50 % reduction, but usually less...
- Risk for re-takes
- Interfering with the AEC
- *Gonad shielding of females can be problematic. Using the other described methods will usually give a higher patient dose reduction*



Conclusion

- Proper selection of equipment is essential
- Good radiographic knowledge, technique and training are crucial for the ALARA approach
- Education and training is important in radiation protection

Thank you!

