

Presented at the NSFS conference, Aug. 26, 2015 (Roskilde, Denmark)



New Danish research laboratory for medical dosimetry

Claus E. Andersen
Jakob Helt-Hansen
Lars R. Lindvold
Bent Lauritzen + 4 PhD's

Center for Nuclear Technologies
Technical University of Denmark
Risø Campus, Build. 201
Frederiksborgvej 399
4000 Roskilde, Denmark

clan@dtu.dk

Duration: 20 min



Outline

Laboratory for Fundamental Medical Dosimetry

1. Radiation sources and research objectives.
2. Specific details (mainly related to calibration system).
3. Radiation protection measurement in pulsed beams.

Donation (21. juni 2012) from:

The John and Birthe Meyer Foundation

Accelerator, cobalt irradiation facility and advanced instrumentation.

The overall objective for the laboratory:

To support development and introduction of new forms of radiotherapy on Danish hospitals and to contribute to the maintenance of high standards for patient safety during modern radiotherapy.

Varian Truebeam linear accelerator system



Full imaging (kV, MV, gating)
Eclipse TPS
Developer mode


Photon beams: 4, 6, 10, 15 and 18 MV
FFF photon beams: 6 and 10 MV
Electron beams: 6, 9, 12, 15, 16, 18, and 20 MeV


Terabalt (UJP Praha, Czech Republic) Cobalt-60 source

455 TBq activity


Fixed horizontal beam

Full computer control







Measurement infrastructure



Other irradiation facilities at DTU Nutech






Types of research at the new laboratory:

1. Detector oriented research


Organic scintillator dosimetry

Alanine dosimetry

3D optical dosimetry



2. Clinically oriented projects



Lung tumor dosimetry



3. Calibrations and standards

Accredited calibrations and measurements

Standards for accelerator dosimetry

Audit methods

Project examples:

Research projects within EURAMET

Small-field dosimetry (Danish Research Council)

Radiobiology

Mikael Jensen

Torsten Groesser



Pil Fredericia

(Hevesy Laboratory)



Medical Radiation Dosimetry Course

www.mrdc.dtu.dk




$$D = \int_{T_{\text{cut}}}^{\infty} \Phi_T \frac{S_{\text{col}}(T)}{\rho} k_{\text{col}}(T) dT$$



The DTU Nutech course "Medical Radiation Dosimetry" provides the theoretical and metrological background for measurement of radiation doses in radiotherapy, diagnostic radiology, and nuclear medicine. The five-day intensive course is based on the classic textbook by F.H. Attix and supplementary material.


Platoon 6 (June 15-19, 2015)



Details about

the cobalt source calibration system

Absorbed dose to water calibration procedure



$$d = 5 \text{ g/cm}^2$$


$$\text{SAD} = 100 \text{ cm}$$

$$A_0 = 10 \times 10 \text{ cm}^2$$

Need cobalt to:

- calibrate one chamber against another (substitution method)
- calibrate chamber against standard

Cobalt calibration rig







Alignment telescope

Characterizing the cobalt source (one example)


Collimator reproducibility





Order of collimator settings (first 4 h)




Order of collimator settings







Sensitivity coefficient: $(0.14 \pm 0.01) \% \text{ pr. mm}$



Radiation protection

in pulsed beams

Motivations:

- (1) NSFS
- (2) previous experience
- (3) no specification from manufacturer






Table 1 Measurement results and standard uncertainties for the two spherical ionization chambers (1 L and 10 L) and the Ludlum 9DP* survey meter in scattered beams from the medical accelerator and the cobalt source at the reference position shown in Fig. 1 and 4. The phantom dose was measured with a thimble chamber.

Beam	Phantom dose Gy
4 MV	4.0
6 MV	8.5
10 MV	9.4
15 MV	9.8
18 MV	10.2
6 MV, FFF	8.0
10 MV, FFF	9.0
Cobalt-60	-
FFF = Flattening Filter Free	

Table 1 Measurement results and standard uncertainties for the two spherical ionization chambers (1 L and 10 L) and the Ludlum 9DP* survey meter in scattered beams from the medical accelerator and the cobalt source at the reference position shown in Fig. 1 and 4. The phantom dose was measured with a thimble chamber.

Beam	Phantom dose Gy	Irradiation time s	Pulse width μ s	Pulse-to-pulse ms	No. of pulses
4 MV	4.0	120	3.7	2.5	48.1×10^3
6 MV	8.5	120	3.7	3.3	36.0×10^3
10 MV	9.4	120	3.5	3.3	36.0×10^3
15 MV	9.8	120	3.2	6.7	18.0×10^3
18 MV	10.2	120	2.9	6.7	18.0×10^3
6 MV, FFF	8.0	150	4.4	4.9	30.9×10^3
10 MV, FFF	9.0	150	4.1	8.3	18.0×10^3
Cobalt-60	-	-	-	-	-

FFF = Flattening Filter Free



Time structure measured
with ME40 + oscilloscope

Table 1 Measurement results and standard uncertainties for the two spherical ionization chambers (1 L and 10 L) and the Ludlum 9DP* survey meter in scattered beams from the medical accelerator and the cobalt source at the reference position shown in Fig. 1 and 4. The phantom dose was measured with a thimble chamber.

Beam	Phantom dose Gy	Irradiation time s	Pulse width μ s	Pulse-to-pulse ms	No. of pulses	H*(10) 1 L μ Sv/h
4 MV	4.0	120	3.7	2.5	48.1x10 ³	260 $\pm$ 30
6 MV	8.5	120	3.7	3.3	36.0x10 ³	480 $\pm$ 50
10 MV	9.4	120	3.5	3.3	36.0x10 ³	380 $\pm$ 40
15 MV	9.8	120	3.2	6.7	18.0x10 ³	910 $\pm$ 100
18 MV	10.2	120	2.9	6.7	18.0x10 ³	1600 $\pm$ 180
6 MV, FFF	8.0	150	4.4	4.9	30.9x10 ³	240 $\pm$ 30
10 MV, FFF	9.0	150	4.1	8.3	18.0x10 ³	180 $\pm$ 20
Cobalt-60	-	-	-	-	-	67 $\pm$ 7

FFF = Flattening Filter Free

Spherical chamber calibrated in Co-60 beam for air kerma, free in air converted to ambient dose equivalent H*(10) using a conversion coefficient of 1.2 Sv/Gy (reservation: energy spectrum and angular dependency).

Table 1 Measurement results and standard uncertainties for the two spherical ionization chambers (1 L and 10 L) and the Ludlum 9DP* survey meter in scattered beams from the medical accelerator and the cobalt source at the reference position shown in Fig. 1 and 4. The phantom dose was measured with a thimble chamber.

Beam	Phantom dose Gy	Irradiation time s	Pulse width μ s	Pulse-to-pulse ms	No. of pulses	H*(10) 1 L μ Sv/h	H*(10) 10 L μ Sv/h
4 MV	4.0	120	3.7	2.5	48.1x10 ³	260 $\pm$ 30	260 $\pm$ 30
6 MV	8.5	120	3.7	3.3	36.0x10 ³	480 $\pm$ 50	500 $\pm$ 60
10 MV	9.4	120	3.5	3.3	36.0x10 ³	380 $\pm$ 40	400 $\pm$ 40
15 MV	9.8	120	3.2	6.7	18.0x10 ³	910 $\pm$ 100	930 $\pm$ 100
18 MV	10.2	120	2.9	6.7	18.0x10 ³	1600 $\pm$ 180	1700 $\pm$ 180
6 MV, FFF	8.0	150	4.4	4.9	30.9x10 ³	240 $\pm$ 30	250 $\pm$ 30
10 MV, FFF	9.0	150	4.1	8.3	18.0x10 ³	180 $\pm$ 20	190 $\pm$ 20
Cobalt-60	-	-	-	-	-	67 $\pm$ 7	67 $\pm$ 7

FFF = Flattening Filter Free



Table 1 Measurement results and standard uncertainties for the two spherical ionization chambers (1 L and 10 L) and the Ludlum 9DP* survey meter in scattered beams from the medical accelerator and the cobalt source at the reference position shown in Fig. 1 and 4. The phantom dose was measured with a thimble chamber.

Beam	Phantom dose Gy	Irradiation time s	Pulse width μ s	Pulse-to-pulse ms	No. of pulses	H*(10) 1 L μ Sv/h	H*(10) 10 L μ Sv/h	H*(10) 9DP* μ Sv/h
4 MV	4.0	120	3.7	2.5	48.1x10 ³	260 $\pm$ 30	260 $\pm$ 30	240 $\pm$ 30
6 MV	8.5	120	3.7	3.3	36.0x10 ³	480 $\pm$ 50	500 $\pm$ 60	440 $\pm$ 40
10 MV	9.4	120	3.5	3.3	36.0x10 ³	380 $\pm$ 40	400 $\pm$ 40	330 $\pm$ 30
15 MV	9.8	120	3.2	6.7	18.0x10 ³	910 $\pm$ 100	930 $\pm$ 100	770 $\pm$ 80
18 MV	10.2	120	2.9	6.7	18.0x10 ³	1600 $\pm$ 180	1700 $\pm$ 180	1500 $\pm$ 150
6 MV, FFF	8.0	150	4.4	4.9	30.9x10 ³	240 $\pm$ 30	250 $\pm$ 30	220 $\pm$ 20
10 MV, FFF	9.0	150	4.1	8.3	18.0x10 ³	180 $\pm$ 20	190 $\pm$ 20	170 $\pm$ 20
Cobalt-60	-	-	-	-	-	67 $\pm$ 7	67 $\pm$ 7	62 $\pm$ 6

FFF = Flattening Filter Free




Ludlum 9DP*

Radiation protection

Dose rate around the accelerator after
4000 MU of 18 MV (40 Gy)


A somewhat unexpected result




Datalogging
1 record pr. 5 s

Positioned after
beam stop


(i.e. instruments were not activated).




Dose rate around the accelerator after 4000 MU of 18 MV (initial 1.5 h)



Dose rate around the accelerator after 4000 MU of 18 MV (after 5 h)



Dose rate around the accelerator after 4000 MU of 18 MV (log scale)



Conclusions



- A laboratory for medical dosimetry has been established. Only few countries have similar facilities.
- Three types of research:
 - (1) Detector projects, (2) Clinical projects, and (3) Calibration and standards (+ radiobiology).
- Showed: Automated cobalt calibration details
- Designed simple test stand for radiation prot. measurements in pulsed beams.
- Found relatively high dose rates due to induced activity after long 18 MV irradiations.